


Boston Chair Dimensions

Please note that depending on the seat cushion built into the chair the seat height may vary slightly from that shown opposite.

Important note – All measurements are approximated as each chair may be slightly different when constructed.



Letter	Feature	Small	Medium	Large
A	Seat Width	457mm (18")	508mm (20")	559mm (22")
B	Overall Width	724mm (28.5")	775mm (30.5")	826mm (32.5")
C	Overall Height	1270mm (50")	1270mm (50")	1270mm (50")
D	Back Height	762mm (30")	762mm (30")	762mm (30")
E	Handle Bar Height	1118mm (44")	1118mm (44")	1118mm (44")
F	Seat Depth	508mm (20")	508mm (20")	508mm (20")
G	Overall Depth Upright	927mm (36.5")	927mm (36.5")	927mm (36.5")
H	Foot Plate Depth	330mm (13")	330mm (13")	330mm (13")
I	Seat Height (Leg Rest to Seat)	394mm (15.5") – 546mm (21.5")	394mm (15.5") – 546mm (21.5")	394mm (15.5") – 546mm (21.5")
J	Arm Height	178mm (7")	178mm (7")	178mm (7")
K	Seat Height (Floor to Seat)	635mm (25")	635mm (25")	635mm (25")
L	Overall Depth Fully Reclined	1829mm (72")	1829mm (72")	1829mm (72")
i	Back Angle	92° – 154°	92° – 154°	92° – 154°
ii	Leg Rest Angle	86° - 170°	86° - 170°	86° - 170°
TIS Range 0°- 16°				